



60 Walnut Ave
Clark, NJ 07066
Tel: 732-882-0200
Fax: 732-882-0300
info@dbtech.com
www.dbtech.com

RAS Build 332 Release Notes

The new build was released in the fourth quarter of 2010.

The new RAS build includes the following new features and enhancements:

RAS*i* General

- Document Types and Document Sets can now be associated with multiple Patient Types. Also, it is no longer required that Document Types and Document Sets be associated with a Patient Type, thus allowing for “global” DocTypes and DocSets.

RAS*i* Viewer

- The RAS*i* Viewer now includes an option to export selected patient documents. The documents are exported to a user-specified folder in PDF format along with a cover/search page to make it easy to navigate and view the exported documents.
- The viewer now supports an application programming interface (API) that allows users to initiate a RAS*i* Patient Folder search from another application.
- The Document Summary pane can now optionally group documents by category.
- The viewer now includes an option to print the contents of the current patient folder or document list. The user can select the data fields that will appear in the printout.
- Patient photos can now be deleted.

RAS*i* Web Viewer

- The RAS*i* Web Viewer now includes scanning support, including Barcode and Batch Scanning.
- The viewer now includes the ability to print Document Sets on demand.
- A Patient Folder Search interface is now included in the viewer.

RAS Report Viewer

- The Report Viewer (including the Web version) can now be configured to support .DOC, .PPT, .RTF, and any other files for which a Windows association exists.
- The Report Viewer (non-Web version) now includes an option to print the contents of the current report list. The user can select the data fields that will appear in the printout.
- The Web Report Viewer can now accept a URL that will take the user directly to a specified report without needing to navigate the report list.

RAS Forms Viewer

- The Forms Viewer now includes an option to print the contents of the current Form’s data list. The user can select the data fields that will appear in the printout.

Management Reports

- A problem where long folder names were truncated on the Folder, Title, and User lists was corrected.